

MAXILLARY MOLAR DISTALIZATION

Patients prefer fixed appliance treatment with no headgear wear!

DISTAL JET™

MOLAR DISTALIZING SYSTEM



A FIXED LINGUAL APPLIANCE THAT UTILIZES NANCE PALATAL ACRYLIC FOR ANCHORAGE AND NiTi COIL SPRINGS WITH AN ADJUSTABLE ACTIVATION LOCK FOR RAPID DISTALIZATION OF UPPER FIRST AND SECOND MOLARS.



DISTAL JET

Appliance offers predictable Unilateral or Bilateral Class II correction without tipping of molars or loss of anterior anchorage, and is easily converted into a passive nance appliance for retention.

- ▶ **SIMPLE** – easy to insert and activate.
- ▶ **ESTHETIC** – the lingual appliance is not visible and is well tolerated by patient.
- ▶ **COOPERATION** – the fixed appliance does not require patient compliance for activation.
- ▶ **RAPID** – NiTi distalizing springs average .91mm per month for short treatment time.
- ▶ **VERSATILE** – a wide variety of design applications available.



Ask for
DISTAL JET
appliances,
fabricated by:

Your full service dental laboratory!
PROTEC
DENTAL LABORATORIES LTD.



34 East 2nd Avenue
Vancouver, BC
Canada V5T 1B1
Main: 604-873-8000
Dispatch: 604-873-1212
Fax: 604-873-8527
Toll Free: 1-800-663-5488

www.protecdental.com info@protecdental.com