



Care and Cleaning of Your Valplast Denture



Ask for Valplast Dentures
fabricated by:

Your full service dental laboratory!

PROTEC
DENTAL LABORATORIES LTD.



Care and Maintenance

Valplast dentures can deteriorate or become lost over time. Therefore, several precautions should be taken to maintain your denture.



**DO NOT SOAK
IN BLEACH**



**DO NOT APPLY
BOILING WATER**



**KEEP IN WATER
WHEN NOT USING**

Wear your denture as instructed. Never attempt to clean your denture by boiling it in hot water. Do not wrap your denture in a napkin or tissue as it could be mistakenly thrown out. Be careful not to leave your denture where it might be stepped on or chewed by a pet.

Regularly check the condition of your denture. Irregularities can irritate teeth and mouth tissues. Be sure to keep your appointments and visit your dentist regularly to have your denture adjusted and your oral health assessed.



Cleaning and Storage

CLEAN ALL SURFACES WITH A SOFT TOOTHBRUSH AND WASH WITH WATER



SOAK DAILY IN VAL-CLEAN OR AN EFFERVESCENT DENTURE CLEANING SOLUTION



STORE IN WATER USING A FIRM, PERFORATED CONTAINER

