




Exclusively
in Canada
from Protec!

iNMAN
aligner
CERTIFIED LAB

Precision Alignment Innovative



iNMANTM aligner

SPRING ALIGNMENT APPLIANCE

Inman Aligner treatment utilizes a corrected set-up arrangement of the anterior teeth and a removable appliance design. Niti coil springs with active lingual and labial components apply gentle steady forces to correct most anterior rotations, crowding, and relapse.

FAST - light forces can rapidly correct alignment in just 6 - 16 weeks

VERSATILE - correct moderate anterior rotations, crowding, and relapse

CONTROL - tooth movement is pre-programmed to treatment set-up

COMFORTABLE - easily inserted and removed as a retainer type appliance



Before



After



PROTEC
DENTAL LABORATORIES LTD.



34 East 2nd Avenue
Vancouver, BC
Canada, V5T 1B1
Tel: 604.873.8000
Toll Free: 800.663.5488
Fax: 604.873.8527

Crown & Bridge • Dentures • Orthodontics
www.protecdental.com info@protecdental.com

