

Mould Number 33° Posteriors



Mould Number 33° Posteriors



Mould Number 0° Posteriors



LS indicates Long Buccal-Short Lingual moulds.

A Lucitone® Denture Resin For Every Need

Lucitone 199®

- Over 40 years of clinical success!
- Outstanding esthetics
- Superior impact resistance
- Premium product = Premium results
- Superb flexural strength
- Can be compression packed or injected
- Microwavable if injected



Lucitone® FRS®



- Nylon flexible partial
- Adjustable flexibility
- Strong, high impact resistance
- Color stability
- Low water absorption
- Can be injected

Lucitone® FasPor™ +

- Available in Lucitone 199 shades
- Color stable polymerization
- Flows easily
- Cures in a pressure pot
- Can be poured or injected



Hy-Pro Lucitone®

- Quick set-up time – ready to pack in 3 minutes
- Up to 60 minutes working time
- Lucitone color and dimensional stability
- Production friendly
- Can be compression packed or injected



Lucitone® Clear

- Excellent physical strength properties
- Nearly equivalent to Lucitone 199 impact strength standards
- Can be compression packed or injected



Characterized Lucitone®

- Highly characterized
- Fibers achieve realistic vascular simulation
- Gingival tone powders available
- Can be compression packed



Lucitone 199® Repair Material



- Excellent bond strength
- Available in all 4 Lucitone 199 shades
- Outstanding overall esthetics
- Can be compression packed

LuciSof® Soft Denture Liner

- Advanced silicone chemistry that stays soft and flexible
- Premium soft denture liner choice
- Long lasting, soft relines
- Available in pre-mixed sheets or injectable cartridges
- Can be compression packed or injected



DENTSPLY

Trubyte®

DENTSPLY International Inc., 570 West College Avenue, York, PA 17405-0872 1-800-786-0085
prosthetics.dentsply.com
Distributed By: DENTSPLY Canada, 161 Vinyl Court, Woodbridge, ON L4L 4A3 1-800-263-1437

©2013 DENTSPLY International Inc. All right reserved.

89564 Rev. 14 (07/13)

Rx only CE

Classic®/New Hue® Mould Chart



DENTSPLY

Trubyte®
Classic®/New Hue®
Plastic Teeth

Trubyte® Classic® & New Hue®

Upper Anteriors

Mould Number	Articulation with lower mould	Width of Central (mm)	Length of Central Without Collar (mm)	Width of 6 Anteriors on Curve/Flat (mm)
133 (2C)		7.65	8.75	45.0 42.0
134 (3D)		8.15	9.25	48.5 44.2
135 (46)		8.87	10.0	52.5 48.0
136 (46)		9.25	10.25	54.5 50.6
2D (2D)		7.85	9.5	48.0 43.0
3D (3D)		8.20	9.25	49.0 44.1
1N (46)		8.40	9.75	49.0 44.7
2N (2N)		7.70	10.25	47.0 42.5
2P (2P)		8.65	11.25	52.0 46.7
A24 (3N)		8.15	9.5	49.0 43.1
A25 (2E)		8.70	9.5	49.5 45.3

Mould Number	Articulation with lower mould	Width of Central (mm)	Length of Central Without Collar (mm)	Width of 6 Anteriors on Curve/Flat (mm)
A26 (46)		9.20	9.75	52.5 47.2
3M (3M)		7.22	9.00	43.0 39.9
3N (3N)		8.25	10.0	48.5 45.0
3P (3P)		8.70	10.5	52.5 48.5
4M (3N)		7.65	9.5	46.0 40.0
4N (2N)		8.40	10.0	49.0 43.9
263 (2N)		7.65	10.0	46.0 42.2
264 (2E)		8.30	10.75	50.0 45.1
266 (26)		9.20	12.0	55.0 49.9
267* (3R)		9.65	13.0	58.0 52.6
268* (3R)		10.20	13.25	60.5 55.1

*Only offered in Classic tooth line.

Lower Anteriors

Mould Number	Articulation with upper mould	Width of Central (mm)	Length of Central Without Collar (mm)	Width of 6 Anteriors on Curve/Flat (mm)
2C (133)		4.60	8.25	33.5 31.9
3M (3M)		4.85	8.25	33.5 31.9
2D (2D)		5.05	9.25	34.0 32.5
2N (2N)		5.12	9.0	35.0 33.1
3D (3D)		5.15	8.75	36.5 34.0
3N (3N)		5.15	9.5	36.5 34.6
2E (A25)		5.60	9.5	37.5 35.5
2P (28)		5.55	9.75	38.5 35.9
46 (A26)		5.65	8.75	39.0 36.7
3P (3P)		5.85	10.0	39.5 38.1

Note: Teeth are shown actual size.

Posteriors

Mould Number	Articulation with upper mould	Width of Central (mm)	Length of Central Without Collar (mm)	Width of 6 Anteriors on Curve/Flat (mm)
26 (266)		5.65	11.5	39.5 38.1
1H (267)		5.90	10.75	41.5 38.1
3R (267)		5.90	11.75	42.5 40.5

Suggested Articulations & Mould Combinations

Upper Mould Number	Articulates with Lower Mould	Combines with 33° Posterior Mould	Combines with 20° Posterior Mould	Combines with 10° Posterior Mould	Combines with 0° Posterior Mould
133	2C	30M	29M	F30	29M
134	3D	30M	31M	F33	31M
135	46	32M	33M	F33	33M
136	46	34M	33M	F33	33M
2D	2D	30L	29L	F32	29M
3D	3D	30M	31M	F33	31M
1N	46	30L	31M	F32	31M
2N	2N	30L	29L	F30	29M
2P	2P	30L	31L	F33	31M
A24	3N	30M	29M	F30	29M
A25	2E	30M	31M	F33	31M
A26	46	30M	31M	F33	31M
3M	3M	30M	29M	F30	29M
3N	3N	30L	31M	F32	31M
3P	3P	34M	33M	F33	33M
4M	3M	30M	29M	F30	29M
4N	2N	30L	31M	F30	31M
263	2N	30M	29L	F32	29M
264	2E	30M	31L	F33	31M
266	26	34L	33M	F33	33M
267	3R	34L	33M	F33	33M
268	3R	34L	31L	F33	33M

Please Note: The 20° mould form dimensions are the only moulds that differ in the New Hue and Classic tooth lines. The Classic moulds are illustrated for you in the pictures above. Below are the New Hue 20° dimensions.

Mould	X	Y	U	L
29M	29.0	32.9	7.7	9.2
29L	29.0	33.2	8.7	9.1
31M	31.0	35.2	8.1	9.9
31L	31.0	35.5	9.0	9.8
33M	33.0	36.9	8.9	10.2

20° Posteriors

29M X=29.0 Y=32.9 U= 8.5 L= 9.3	
29L X=29.0 Y=33.1 U= 8.4 L= 9.4	
31M X=31.0 Y=35.1 U= 9.1 L=10.0	
31L X=31.0 Y=35.5 U= 8.9 L= 9.8	
33M X=33.0 Y=37.4 U= 9.3 L=10.4	


10° Posteriors

F30 X=30.0 Y=34.0 U= 8.6 L= 8.3	
F32 X=32.0 Y=33.6 U= 8.7 L= 9.3	
F33 X=33.0 Y=35.0 U= 9.5 L= 9.7	

Shade Availability

New Hue Anterior and Posterior Shades:
59, 61, 62, 65, 66, 67, 69, 77, 81, 87
(New Hue 10-tab Shade Guide Available)

Classic Anterior and Posterior Shades*:
CW2™, A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4, 59C, 62C, 65C, 66C, 67C, 69C, 77C, 81C (Classic 25-tab Shade Guide Available)



The average width of a 1 x 4 set of upper posteriors = X;
the lower posteriors = Y.
The average depth of the upper left first molar = U;
the lower left first molar = L.
Processing instructions are available upon request.

*The A1-D4 designations correspond to the Vita® Classical Shade Guide which is meant to be a guide, not a match.
Vita is not a registered trademark of DENTSPLY International Inc.

ISO22112: 2005
ANSI/ADA Specification no. 15-2008